

# TITRIVIN

## A range of wine standards to meet your needs

### ♦ The calibration reference :

⇒ **TITRIVIN AA** series supplies a range of ordinate values for the 10 most common parameters to calibrate the response patterns of automatic equipment

⇒ **TITRIVIN IR** series is a calibration range (ascending alcohol, descending sugar) designed to measure the reducing sugars/alcohols pair by infrared reflectance in fermenting wines.

*1 box of these standards is 20 x 10 mL ampoules*

### ♦ The control reference :

Up to 16 parameters

⇒ **TITRIVIN BTA** used every series samples of wine analysed. If a deviation in value measured is observed over time then remedial action is called for

*1 box = 20 x 10 mL ampoules*

⇒ **TITRIVIN BTB** is designed for checking manual analyses

*1 box = 6 x 240 mL bottles*

Find reference values on our website

[www.titrivin.com](http://www.titrivin.com)

FRANCE



VINE AND WINE DEPARTMENT  
39 rue Michel Montaigne – BP115  
F-33294 BLANQUEFORT CEDEX  
Tél. + 33 (0) 5 56 35 00 00  
Fax + 33 (0) 5 56 35 58 59

[titrivin@titrivin.com](mailto:titrivin@titrivin.com)


[www.titrivin.com](http://www.titrivin.com)

FRANCE AND ABOARD

Laboratoires Dujardin-Salleron  
18 rue Henri Barbusse  
F-94117 ARCUEIL CEDEX  
Tél. + 33 (0) 1 45 46 04 05  
Fax + 33 (0) 1 15 46 01 13

[info@dujardin-salleron.com](mailto:info@dujardin-salleron.com)

[www.dujardin-salleron.com](http://www.dujardin-salleron.com)



## User's guide

# TITRIVIN reference materials



## In your laboratory, TITRIVIN reference materials are essential....

### ♦ To calibrate equipment :

On a certificate of analysis, you find, for a parameter, a **reference value ( $V_{ref}$ ) with its uncertainty ( $u$ )**.

This means the "true" value is between  $[V_{ref} - u ; V_{ref} + u]$  called the **confidence interval**.

Most of the time, you will use the reference value to calibrate.

But **to improve your calibration curve**, you could use **a value within to the confidence interval**.



### ♦ To carry out quality control :

In laboratories, there are two types of quality control :

#### ➤ External quality control :

The laboratory takes part in inter-laboratory analytical comparisons organised by for example the Oenologues de France or BIPEA...

#### ➤ Internal quality control :

The laboratory monitors the quality of their results using a control chart. This control is easily carried out with TITRIVIN wine standards.

When you use TITRIVIN to monitor your control chart, **your target value**, that is to say the average of your results **must be within the confidence interval**.



### ♦ To assess uncertainties

This is a direct use of the quality control with TITRIVIN

When your control chart is in place, you can calculate your reproducibility standard deviation :

$$\sigma_{repro}$$

Therefore, you can assess your uncertainty from this range of analytical results.

Uncertainties can be evaluated by :

$$\text{Uncertainty} = 2 \times \sqrt{\sigma_{repro}^2 + \frac{u^2}{2}}$$